Dansk	2
Deutsch	10
English	19
Español	27
Ελληνικά	36
Français	45
Italiano	54
Nederlands	63
Norsk	72
Polski	80
Português	89
Svenska	98

Velkommen

Tillykke med din nye Osstell® IDx.

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Tel.: +46 31 340 82 50
 Fax: +46 31 413 115
 E-Mail: support@osstell.com

Welcome

Congratulations on the purchase of your new Osstell® IDx.
Before you start using your Osstell IDx, please read through the entire manual.

Table of content

1) Table of content	19
2) Warnings and Precautions	20
3) Indications for use	20
4) Description	20
5) Safety Symbols	21
6) Before You Start	21
7) Preparation and Set-up	21
8) How to Measure	22
9) View Measurement	23
10) Measure in the Classic Mode	23
11) View Measurement in Classic Mode	24
12) Interpret the Result	24
13) Lock the unit	24
14) Cleaning and Maintenance	24
15) Firmware update	25
16) Backup	25
17) Technical information	25
18) Troubleshooting	26
19) Support	26

2) Warnings and Precautions

Warning: The Probe emits an alternating magnetic field that potentially could interfere with cardiac pacemakers!

Precautions:

- Do not use the device in the presence of explosive or flammable materials.
- To avoid interference with other equipment, the probe should not be held close to electronic devices.
- Use only the power supply provided with the instrument.
- SmartPegs™ are for single use only.
- The Osstell IDx Probe and the SmartPeg Mount must be sterilized before use.

3) Indications for use

Osstell IDx is indicated for use in measuring the stability of implants in the oral cavity and maxillofacial region. Intended user profile: clinical professionals. No special experience/training is needed.

4) Description

The Osstell IDx has a built-in patient database and enables the user to connect measurements to implant position and patient through the touch screen. The Osstell IDx involves the use of the noninvasive technique, Resonance Frequency Analysis for measuring implant stability in the oral cavity and the maxillofacial region. The Osstell IDx system includes the use of a SmartPeg which is attached to an implant or abutment using the SmartPeg Mount. The SmartPeg is activated by magnetic pulses from the probe, which is connected to the instrument. The resonance frequency, which is the measure of implant stability, is calculated from the response signal. Results are presented on the instrument display as an ISQ value (Implant Stability Quotient), which is scaled from 1 to 100. The higher the number, the greater the stability. The first measurement is taken at implant installation to get a baseline for future measurements. Before the final restoration, a second measurement is taken which makes it possible to observe the stability of the implant. With the graphical view, the implant stability development can be monitored. The software can be updated online thanks to internet connectivity. Osstell IDx can add important information to the evaluation of implant stability and can be used as part of an overall treatment evaluation program. The final implant treatment decisions are the responsibility of the clinician.

The SmartPeg is available with different connection geometries to fit all major implant products on the market. On osstell.com you can find all available SmartPegs and in the Osstell app.

Your Osstell IDx system includes the following items:

- ① Osstell IDx Instrument
- ② Osstell ISQ Probe
- ③ Osstell SmartPeg Mount
- ④ Osstell TestPeg
- ⑤ Osstell Charger
 - Osstell Plug Adapter
 - Osstell IDx Manual
 - Osstell USB

Applied Parts:
Osstell ISQ probe

5) Safety Symbols

	Follow instructions for use
	Class II equipment
	Type BF equipment
IP20	Protected against solid foreign objects of 12.5mm \varnothing and greater. No protection against water.
	Do not reuse
LOT	Lot/batch number
SN	Serial number
STERILE R	Sterile / Sterilization by radiation
	Manufacturer
	Caution
	Use-by Date
	Sterilizable up to 135 degrees Celsius
	Li-ion battery
RxOnly	Prescription use only. U.S. Federal Law restricts this device to sale by or on the order of a licensed practitioner. - <i>Is not intended for the Canadian market.</i>

6) Before You Start

Battery

Connect the power supply to the instrument at the inlet marked with .

Use the adapter to the power supply that is compatible with your wall sockets.

The rechargeable battery should be charged for at least 3 hours before its first use.

A blue light on the on/off key indicates charging when the instrument is off.

Register your Device

You must register your Osstell IDx. The serial number can be found on the back of the instrument, as well as in the settings section of the software.

For registration assistance, please visit osstell.com/support.

7) Preparation and Set-up

Start the instrument by pressing the on/off key .

The instrument is turned off by pressing the on/off key.

Enter the settings menu by pressing .

Settings

Audio	Sound on/off: Activate by touching the circle. Volume: To increase volume, press +. To decrease volume, press -
Date and Time	Time: Enter the time, press OK. Date: Enter the year, month and day, press OK. Time zone: Choose in scroll bar.
Security	Change of pincode: Press ****, Enter 4-digit code, OK.
Language	Change of language: Choose in scroll bar.

Power

Display	To increase brightness, press + To decrease brightness, press -
Save power	Time to standby: Enter the time when the instrument is to go into standby mode when not in use, and end with OK. To activate - touch the screen. Time to off: Enter the time when the instrument is to shut down after not being used, and end with OK.

Exam

Work mode	Touch the circle and Classic mode is selected.
Dental notation	Choose UNS or FDI dental notation mode using the arrows.
Reset Classic	To delete all measured values in Classic mode: Press OK.

Information

Versions	Version designation of the software
Hardware	Serial number. Seal: Reset time, date, time zone, language, pincode
Register product	Unlock: Enter PUK code, press OK.

Net service

Connection	Network enabled: Activate by touching the circle. WiFi name: Select your wifi network and enter your wifi password.
Net backup	Auto backup: Activate by touching the circle. Manage backup: Export and import of data. New backup password: Press OK. Enter new password - press OK.

TestPeg

The TestPeg may be used for testing and learning how to use the system. Place the TestPeg on the table or hold it in your hand. Turn on the instrument and hold the probe approximately 3-5mm away from the top of the TestPeg.

8) How to Measure

1. Connect the probe to the instrument.
2. Start the instrument .
3. Enter patient's mode .
4. Select New patient to create a new patient data card, type in the required information and end with Save.
5. Select the position where the implant has been inserted by touching the correct tooth position. Dental notation can be changed, see settings.
6. Tie a piece of string or thread through the SmartPeg Mount. The thread should be placed around the wrist when the SmartPeg Mount is used to avoid dropping it.
7. Place the SmartPeg into the SmartPeg Mount  ①. The SmartPeg is magnetic, and the SmartPeg Mount will hold the SmartPeg as it is carried to the implant. Attach the SmartPeg to the implant or abutment by screwing the SmartPeg Mount using finger force of approximately 4-6 Ncm  ②. Do not over-tighten, to avoid destroying the SmartPeg threads.

8. Hold the probe close (3-5mm) to the top of the SmartPeg without touching it ③. Measure in both the mesiodistal ④ and buccolingual direction ⑤. Do not measure from above ⑥. An audible sound indicates that a measurement has been performed and the value, time and date are shown on the display. Measured values can be deleted by touching the trash can . Note that values cannot be deleted after saving.
9. Remove the SmartPeg by using the SmartPeg Mount.

Note! The SmartPeg is for single use only!
10. Save the data by “End Patient”.
11. The instrument is turned off by pressing the on/off key for 2 seconds .
12. Remove the probe from the instrument by gently pulling the probe connector straight out.

Do not pull the cable.

9) View Measurement

Enter the patient list, choose patient and click on the implant to view. A blue circle indicates the selected implant and the tooth number is shown above. All measured data is found on the display. The first measured ISQ value is shown to the left. To show the remaining ISQ values, slide your finger from right to left. Below each ISQ value, the date of measurement is shown. All values are presented in a graphical view: to show, choose .

10) Measure in the Classic Mode

For how to choose Classic mode, see Exam mode section.

1. Connect the Probe to the instrument.
2. Start the instrument .
3. Tie a piece of string or thread through the SmartPeg Mount. The thread should be placed around the wrist when the SmartPeg Mount is used.
4. Connect a SmartPeg to the SmartPeg Mount ①. The SmartPeg is magnetic, and the SmartPeg Mount will hold the SmartPeg as it is carried to the implant. Attach the SmartPeg to the implant or abutment by screwing. Use approximately 4-6Ncm of torque ②. Do not over-tighten, to avoid destroying the SmartPeg threads.
5. Hold the probe close (3-5mm) to the top of the SmartPeg without touching it ③. Measure in both the mesiodistal ④ and buccolingual direction ⑤. Do not measure from above ⑥. An audible sound indicates that a measurement has been performed and the value, time and date are shown on the display.
6. Remove the SmartPeg using the SmartPeg Mount. **Note! The SmartPeg is for single use only!**
7. The instrument is turned off by pressing the on/off key for 2 seconds .
8. Remove the probe from the instrument by gently pulling the probe connector straight out.

Do not pull the cable.

11) View Measurement in Classic mode

The ISQ values are shown directly on the screen. To browse between the measured values, use the arrows.

How to Measure on an Abutment

When a measurement is made on an abutment or on an implant with a “built-in” abutment, the ISQ value will be lower compared to a measurement made on the implant. This is due to the difference in height above bone. To find out the ISQ difference to the measurement performed at implant level, a measurement should be taken on the implant before the abutment is attached and then a second measurement on the abutment.

12) Interpret the Result

Implant Stability

An implant has different degrees of stability in different directions. The total stability consists of the implant stability in relation to the surrounding bone, and the stability of the bone itself. There is always a direction where the stability is lowest, and a direction where the stability is highest and these two directions are perpendicular to each other. The SmartPeg measures the stability in both directions, thus two different ISQ values can be achieved on the same implant. Sometimes the two ISQ values will be very close to each other, or even the same. The high value, which in most cases is found in the mesiodistal direction, mainly reflects the stability in relation to the bone. If a lower value is found, it reflects more of the total stability, where the bone anatomy is a factor.

The ISQ Value

Stability measurements using the Osstell IDx should be performed after the implant or abutment has been placed, assuming there is access to the implant. Measurements should be performed at implant placement and before the implant is loaded or the abutment is connected. After each measurement, the ISQ values are recorded and used as the baseline for the next measurement performed. A change in the ISQ value reflects a change in implant stability.

In general, an increase in ISQ values from one measurement time to the next indicates a progression towards higher implant stability while a decrease in ISQ values indicates a loss in stability and, possibly, implant failure. A stable ISQ value would indicate no change in stability. ISQ values have not been correlated with other methods of implant mobility measurements.

Osstell Connect

Login to [OsstellConnect.com](https://www.osstellconnect.com) for access to tools to interpret and analyze your data. Patient names are not uploaded to Osstell Connect for security reasons, so to find the corresponding patient use “Reference number” as an identifier.

13) Lock the unit

Press the on/off key  to lock the instrument. To unlock the unit: enter the PIN code.

14) Cleaning and Maintenance

Instrument:	can be cleaned using a damp cloth with water or isopropyl alcohol solution.
SmartPegs:	single use only.
Probe and SmartPeg Mount:	must be autoclaved according to the recommended sterilization method, validated to sterility assurance levels (SAL), according to ISO 17665-1 and ISO 17664. Other sterilization cycles may also be suitable, although individual clinicians or hospitals are advised to validate other methods before use. The Probe and SmartPeg Mount should be placed in suitable packaging for the sterilization process.

Sterilization

<p>Method: gravity displacement steam Exposure temperature: 135 °C (275 °F) Exposure time: 10 minutes Drying time: 30 minutes</p>	<p>Warnings: do not exceed 135 °C (275 °F). Remove from autoclave immediately after finished cycle. Allow to cool down before use; handle carefully when hot. Do not wash in dishwasher. Do not pour liquids directly into any of the connectors.</p>
--	--

Waste

The Osstell IDx (instrument, probe, charger and USB cable) should be recycled as electrical equipment. SmartPegs should be recycled as metal. Whenever possible, the battery should be disposed in a discharged state to avoid heat generation through inadvertent short-circuiting.

15) Firmware update

The instrument firmware can be upgraded in two ways:

1. Connect the instrument to the internet, (see SETTINGS section) and an automatic search for new updates will start. If a new version has been released, a pop-up window will appear and instruct you on how to proceed. or
2. Download the latest version on **OsstellConnect.com** to a USB memory. Insert the USB memory into the USB port of the Osstell IDx while it is on. The update will start automatically. Follow the instructions to finish.
Attention! All other equipment (not provided with the Osstell IDx) that is connected to the Osstell IDx must meet the applicable IEC standard.

16) Backup

A back-up of the Osstell IDx data is automatically saved if you are connected to WiFi or a network. It can also be saved to a USB Memory stick. Go to settings and “export” to create a backup. Choose “import” to restore an existing backup.

17) Technical information

Technical description

Osstell IDx is CE-marked according to MDD in Europe. Osstell IDx is constructed according to, and fulfills, the standards EN 60601-1 (Class II, internally powered, type BF applied parts. Not AP or APG equipment. Not protected against ingress of water), and EN 60601-1-2. The symbols used, follow the European standard EN 60601-1 and ISO standards 9687 and 15223 as far as possible.

Battery and charger

The instrument contains a lithium-ion battery.

The instrument should only be charged with the Osstell IDx power supply. Battery charging is indicated by a flashing blue LED. The battery symbol displays the approximate status of the battery. A filled symbol indicates batteries are at 100% capacity.

The power supply operates with mains voltage from 100 to 240 VAC, 50-60 Hz. Thus it can be used in most countries in the world using the appropriate plug adapter. Use only in dry, indoor locations.

To disconnect the equipment from supply mains, disconnect the charger from the mains socket outlet.

Accuracy

ISQ accuracy is within +/- 0.5 ISQ units for a single SmartPeg. The accuracy between different SmartPegs is +/- 2 ISQ units.

Power, weight & size:

Rated power: 12V, type FW 7556M/12
 Instrument size: 203mm x 163mm x 72mm
 Package size: 270mm x 243mm x 105mm
 Instrument weight: 0.85kg
 Gross weight: 2.1kg

Environmental conditions during transport:

Temperature: -40°C to +70°C
 Relative humidity: 10% to 95%
 Pressure: 500 hPa to 1060 hPa

Environmental conditions during use:

Temperature: +10°C to +40°C
 Relative humidity: 30% to 75%
 Pressure: 700 hPa to 1060 hPa
 IP class: IP20

18) Troubleshooting

Problem	Cause	Possible Solution
No measurement	Probe not connected.	Connect the probe.
	Electromagnetic interference.	Remove the source of electromagnetic interference.
	Probe is held too far away from the SmartPeg.	Hold the probe 3-5mm away from the SmartPeg.
Charger does not fit the electrical outlets	Wrong Osstell plug adapter.	Use correct Osstell plug adapter.
Instrument does not start	Uncharged battery.	Charge the Osstell IDx.
Difficult to measure in an exact buccolingual direction	No space.	Try to measure at a slightly different rotational angle.
Difficulties screwing in the SmartPeg	Wrong SmartPeg.	Ensure that the SmartPeg is compatible with the implant system.
The instrument does not work	The Osstell IDx has not been registered.	Register the instrument at OsstellConnect.com or via the IDx. Follow the step by step instructions on osstell.com/support
	The probe is not connected properly.	Connect the probe.
The tooth view is gone	Classic examination mode is activated.	Deactivate classic mode, see SETTINGS section.

If your problem continues, please contact your distributor or our support team for further assistance. In the event of an instrument malfunction, the Osstell IDx and accessories should be sent to the manufacturer for repair.

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Osstell AB