Dansk	2
Deutsch	14
English	26
Español	38
Ελληνικά	50
Français	62
Italiano	74
Nederlands	86
Norsk	98
Polski	110
Português	122
Svenska	134

# Velkommen

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Li-ion

# Welcome

Congratulations on the purchase of your new Osstell® IDx. Before you start using your Osstell IDx, please read through the entire instructions for use.

# Table of content

1)	Introduction	27
2)	Warnings and Precautions	27
3)	Intended Use	27
4)	Indications for Use	27
5)	Description	28
6)	Safety Symbols	28
7)	Before You Start	30
8)	Preparation and Set-up	30
9)	TestPeg	31
10)	SmartPeg™	31
11)	How to Measure	31
12)	View stored Measurements	32
13)	How to Measure on an Abutment	32
14)	How to Measure in Classic Mode	32
15)	View Measurement in Classic Mode	32
16)	Interpret the Result	33
17)	Data Connection to Osstell Connect	33
18)	Lock the unit	33
19)	Cleaning and Maintenance	33
20)	Firmware update	34
21)	Backup	34
22)	Technical information	34
23)	Troubleshooting	36
24)	Service and Support	37
25)	Waste and Disposal	37

# 1) Introduction

#### Qualifications of the user

This medical device is intended to be used by qualified dentists, doctors, surgeons, or specialist staff appointed by the responsible clinician.

#### Responsibilities of the user

Read through the entire instructions for use before using this device.

Observe the Warnings and Precautions.

Any serious incident that has occurred in relation to the medical device should be reported to the manufacturer and the competent authority!

#### Responsibility of the manufacturer

The manufacturer can only accept responsibility for the safety, reliability and performance of the medical device when compliance with the following instructions is ensured:

- The medical device must be used in accordance with these instructions for use.
- Modifications or repairs may only be undertaken by Osstell.
- Unauthorized opening of the unit invalidates all claims under warranty and any other claims.

In addition to unauthorized disassembly, modification or repair of the unit and non-compliance with these instructions for use, improper use will void the warranty and release Osstell from all other claims.

### 2) Warnings and Precautions

#### Warnings

▲ The Probe emits an alternating magnetic field that potentially could interfere with cardiac pacemakers! Keep the probe away from implanted electronic devices. Do not place the probe on the patient's body.

- ⚠ Use only the power supply provided with the instrument.
- ▲ The SmartPegs are disposable and should only be used for one or multiple measurements at one treatment session, for a single patient only (to avoid cross-contamination). Repeated re-use may result in false readings due to wear and tear of the soft aluminium SmartPeg threads. Do not use if the product sterile barrier system or its packaging is compromised.
- ⚠ The Osstell ISQ Probe and the SmartPeg Mount must be sterilized before use.

#### Precautions

A Do not use the device in the presence of explosive or flammable materials.

- $\Lambda$  To avoid interference with other equipment, the probe should not be held close to electronic devices.
- ▲ To remove the probe from the instrument, pull gently holding around the probe connector. Do not pull the cable and do not twist the connector.

### 3) Intended Use

Osstell IDx is intended for use as a Dental Implant Stability Analyzer.

### 4) Indication for Use

Osstell IDx is indicated for use in measuring the stability of implants in the oral cavity and maxillofacial region.

#### Conditions

Surgically placed implants or abutments for which there is space to attach a compatible SmartPeg.

#### Reasons for use

The Osstell IDx can add important information to the evaluation of implant stability and can be used as part of an overall treatment evaluation program. The final implant treatment decisions are the responsibility of the clinician.

#### Contraindications

Osstell IDx is contraindicated for implant systems to which the SmartPeg could not be attached for mechanical incompatibility reasons. See section 10 for more information about SmartPegs. Osstell IDx is contraindicated when used together with Pegs not approved by Osstell. Osstell IDx is contraindicated where it is not possible to attach the SmartPeg due to lack of space, or where it impinges on other artificial or anatomical structures.

### 5) Description

The Osstell IDx is a handheld instrument that involves the use of the non-invasive technique, Resonance Frequency Analysis. The system involves the use of a SmartPeg attached to the dental implant or abutment by means of an integrated screw. The SmartPeg is excited by a magnetic pulse from the measurement probe in the handheld instrument. The resonance frequency, which is the measure of implant stability, is calculated from the response signal. Results are displayed on the instrument as the Implant Stability Quotient (ISQ). The ISQ is scaled from 1 to 100. It is a measurement of the stability of the implant, and is derived from the resonance frequency value obtained from the SmartPeg. The higher the number, the greater the stability.

The software can be updated online by internet connectivity. The Osstell IDx has a built-in patient database and enables the user to connect measurements to implant position and patient through the touch screen.

#### Your Osstell IDx system includes the following items:

- 1 Osstell IDx Instrument
- 2 Osstell ISQ Probe
- ③ Osstell SmartPeg Mount
- ④ Osstell TestPeg
- ⑤ Osstell IDx Power Supply
- (6) Osstell IDx Instructions for Use Osstell Plug Adapter

Applied Parts: Osstell ISQ Probe

# 6) Safety Symbols

8	Follow instructions for use
	Consult instructions for use
	Class II equipment
$\square$	Caution
	See section 2) Warnings and Precautions
Ŕ	Type BF applied part
	Manufacturer
VYYY-MM-DD	Country and date of manufacture
SN	Serial number
Li-ion	Do not dispose of with domestic waste. Li-ion battery.

<b>CE</b> 0297	CE mark with identification number of the notified body.
IP20	Protected against solid foreign objects of 12.5 mm Ø and greater. No protection against water.
1	Temperature limitation
<u>%</u>	Humidity limitation
Í	Atmospheric pressure limitation
135°C 555	Sterilizable up to 135 degrees Celsius
2	Do not reuse
$\square$	Use by date
LOT	Lot/batch code
STERILER	Sterilized using irradiation
$\mathbf{R}_{\mathbf{k}}$ only	For US market only: Prescription use only. U.S. Federal Law restricts this device to sale by or on the order of a licensed practitioner.
تظمر	Only use with power supply type FW8001M/12
V	Supply voltage of the control unit
Α	Supply current
Hz	Frequency of the alternating current
c <b>Al</b> us	UL Component Recognition. Mark indicates compliance with Canadian and U.S. requirements
REF	Catalog Number
MD	Medical device
	Data Matrix code for product information including UDI (Unique Device Identification)
	Do not use if package is damaged and consult instructions for use
Ť	Keep dry
类	Keep away from sunlight
	Single sterile barrier system with protective packaging outside
$\bigcirc$	Single sterile barrier system

### 7) Before You Start

### Battery

Connect the power supply on the back of the instrument. Use the plug adapter to the power supply that is compatible with your wall socket. The rechargeable battery should be charged for at least 3 hours before its first use. A blue light on the on/off key indicates charging when the instrument is off.

### **Register your Device**

You must register your Osstell IDx. The serial number can be found on the back of the instrument, as well as in the settings section of the software.

For registration assistance, please visit **osstell.com/get-started** 

### 8) Preparation and Set-up

Start the instrument by pressing the on/off key **O**. The instrument is turned off by pressing the on/off key for 2 seconds. Enter the settings menu by clicking **O**. Select More to enter the different settings.

Settings

Audio	Sound on/off: Activate by clicking the circle. Click again to de-activate. Volume: To increase volume, click +. To decrease volume, click –
Date and Time	Time: Enter the time, click OK. Date: Enter the year, month and day, click OK. Time zone: Choose in scroll bar.
Security	Change of pincode: Click ****. Enter 4-digit code, OK.
Language	Change of language: Choose in scroll bar.

#### Power

Display	To increase brightness, click + To decrease brightness, click -
Save power	Time to standby: Enter the time when the instrument is to go into standby mode when not in use, and end with OK. To activate - touch the screen. Time to off: Enter the time when the instrument is to shut down after not being used, and end with OK.

#### Exam

Work mode	Click the circle and Classic mode is selected.
Dental notation	Choose UNS or FDI dental notation mode using the arrows.
Reset Classic	To delete all measured values in Classic mode: Click OK.

#### Information

Versions	Version designation of the software and the SmartPeg list.
Hardware	Serial number. Factory Reset: All patient, implant and measurement data will be deleted on the IDx and Osstell Connect. All settings will be set to factory default.
Register product	Unlock: Enter PUK code, press OK.

### Net service

Connection	Network enabled: Activate by clicking the circle. Click again to de-activate. WiFi name: Select your wifi network and enter your wifi password.		
Net backup	Auto backup: Activate by clicking the circle. Click again to de-activate.		
	Manage backup: Export and import of data.		
	New backup password: Press OK. Enter new password - press OK.		

### 9) TestPeg

The Osstell TestPeg may be used for testing and learning how to use the system. Place the Osstell TestPeg on the table or hold it in your hand. Turn on the instrument, select a tooth position and hold the probe tip approximately 3-5 mm away from the top of the Osstell TestPeg. The Osstell IDx should start to measure and present an ISQ value of 55 + /-2 ISQ.

### 10) SmartPeg™

The SmartPeg is available with different connection geometries to fit all major implant products on the market. You can find all available SmartPegs on **osstell.com/smartpegguide**.

The SmartPegs are disposable and should only be used for one or multiple measurements at one treatment session, for a single patient only (to avoid cross-contamination). Repeated re-use may result in false readings due to wear and tear of the soft aluminium SmartPeg threads.

### 11) How to Measure

A first measurement should be taken at implant placement to get a baseline for future measurements throughout the healing process. Before the final restoration, a measurement is taken which makes it possible to observe the stability development of the implant. We recommend you studying the more detailed information (videos and quickguides) available on **osstell.com/get-started**, to utilize the full functionality of your IDx.

- 1. Connect the probe to the instrument.
- 2. Start the instrument **O**.
- 3. Press and hold a tooth to prepare site.
- 4. Enter data about the patient. We recommend using the same patient number as you use in your own patient management system **2**.
- Place the SmartPeg into the SmartPeg Mount ①. The SmartPeg i magnetic, and the SmartPeg Mount will hold the SmartPeg. Attach the SmartPeg to the implant or abutment by screwing the SmartPeg Mount using finger force of approximately 4-6 Ncm ②.
  Do not over-tighten, to avoid damaging the SmartPeg threads.

6. Hold the probe close (3-5 mm) to the top of the SmartPeg without touching it ③. Measure in both the Buccal-Lingual ④ and Mesial-Distal direction ⑤. Do not measure from above ⑥. An audible sound indicates that a measurement has been performed and the value, time and date are shown on the display. Measured values can be deleted by clicking the trash can. Note that values cannot be deleted after saving.

- 7. Save the data by clicking "Save".
- 8. Remove the SmartPeg by using the SmartPeg Mount.
- 9. The instrument is turned off by pressing the on/off key for 2 seconds 🔘.
- 10. To remove the probe from the instrument, pull gently straight out, holding around the probe connector. Do not pull the cable and do not twist the connector.

### 12) View stored Measurements

Enter the patient list **I**, choose patient and click on the implant to view. A blue circle indicates the selected implant and the tooth number is shown above. All measured data is found on the display. The first measured ISQ value is shown to the left. To show the remaining ISQ values, slide your finger from right to left. Below each ISQ value, the date of measurement is shown. The values can be graphed by clicking **I**.

### 13) How to Measure on an Abutment

When a measurement is made on an abutment or on an implant with a "built-in" abutment, the ISQ value will be lower compared to a measurement made on the implant. This is due to the difference in height above bone level. To find out the ISQ difference to the measurement performed at implant level, a measurement should be taken on the implant before the abutment is attached and then a second measurement on the abutment.

### 14) How to Measure in Classic Mode

To choose Classic mode, see section 8) Preparation and Set-up, Exam mode section.

- 1. Connect the Probe to the instrument.
- 2. Start the instrument **O**.
- Place the SmartPeg into the SmartPeg Mount ①. The SmartPeg is magnetic, and the SmartPeg Mount will hold the SmartPeg. Attach the SmartPeg to the implant or abutment by screwing the SmartPeg Mount using finger force of approximately 4-6 Ncm ②.
  Do not over-tighten, to avoid damaging the SmartPeg threads.
- 4. Hold the probe close (3-5 mm) to the top of the SmartPeg without touching it ③. Measure in both the Buccal-Lingual ④ and Mesial-Distal direction ⑤. Do not measure from above ⑥. An audible sound indicates that a measurement has been performed and the value, time and date are shown on the display. Measured values can be deleted by clicking the trash can. Note that values cannot be deleted after saving.
- 5. Remove the SmartPeg using the SmartPeg Mount.
- 6. The instrument is turned off by pressing the on/off key for 2 seconds 🚳
- 7. To remove the probe from the instrument, pull gently straight out, holding around the probe connector. **Do not pull the cable and do not twist the connector.**

### 15) View Measurement in Classic mode

The ISQ values are shown directly on the screen. To browse between the measured values, use the arrows.

1

### 16) Interpret the Result

### Implant Stability

An implant can have different stability in different directions. The stability of the implant is dependent of the surrounding bone configuration. There is often a direction where the stability is lowest, and a direction where the stability is highest and these two directions are most often perpendicular to each other.

To find the lowest stability (lowest ISQ value) it is recommended to measure from two different directions. The lowest stability is in most cases found in the buccal-lingual direction. The highest stability is in most cases found in the mesial-distal direction.

#### The ISQ Value

Assuming there is access to the implant, ISQ measurements should be performed at implant placement and before the implant is loaded or the abutment is connected. After each measurement, the ISQ values are used as the baseline for the next measurement performed. A change in the ISQ value reflects a change in implant stability. In general, an increase in ISQ values from one measurement time to the next indicates a progression towards higher implant stability while a decrease in ISQ values indicates a loss in stability and, possibly, implant failure. A stable ISQ value would indicate no change in stability.

The ISQ scale clinical guidelines are available by clicking on this symbol on the IDx screen. 🗐

### 17) Data Connection to Osstell Connect

Osstell Connect (osstellconnect.com) is an online service for data transfer, storage, display and overview of your data. You can have your Osstell IDx connected to Osstell Connect by WIFI or ethernet cable.

Before you can start using data connection to Osstell Connect you have to register your Osstell IDx. The serial number can be found on the back of the instrument. For registration assistance and features of data connection to Osstell Connect, please visit: **osstell.com/get-started**.

### 18) Lock the unit

Press the on/off key 🛛 to lock the instrument. To unlock the unit: press the on/off key and enter the PIN code.

### 19) Cleaning and Maintenance

Before each use, parts in contact with the patient (Probe, Probe cable and SmartPeg Mount) should be cleaned and sterilized by autoclaving. The instrument does not require regular maintenance. In the event of an instrument malfunction contact the local sales representative or distributor for further instructions.

Instrument:	Can be cleaned using a damp cloth with water or isopropyl alcohol solution. Do not use other cleaning fluids as they can damage the top cover of the instrument.
SmartPegs:	Single use only.
TestPeg:	Is not used intra-orally, does not require sterilization.
Probe and Smart- Peg Mount:	Must be autoclaved according to the recommended sterilization method, validated to sterility assurance levels (SAL), according to ISO 17665-1 and ISO 17664. The Probe and SmartPeg Mount should be placed in suitable packaging for the sterilization process. Loosely coil the cable to allow for packaging without damage to the cable.

### Sterilization

Method: Pre-vacuum	<b>Exposure</b> <b>temperature</b> 132 °C (270 °F)	Exposure time 4 min 2 min	Carefully inspect the probe and cable for damage or wear. Ha wash the probe using a neutral instrument detergent. Do not submerge the cord. Rinse and dry; carefully inspect the probe and cord for damage and wear.
Gravity	134 °C (273 °F)	10 min	Sterilize the probe and cable according to the autoclave manu- facturer's instructions. Do not wash in dishwasher.
<b>Warnings:</b> do not exceed 135 °C (275 °F). <b>Drying time:</b> 30 minutes			Do not pour liquids directly into any of the connectors. Store sterile goods dust-free and dry.

### 20) Firmware update

#### The instrument firmware can be upgraded in two ways:

- Connect the instrument to the internet, see section 8, Preparation and Set-up, Settings section. An automatic search for new updates will start. If a new version has been released, a pop-up window will appear and instruct you on how to proceed. Or;
- 2. Download the latest version on OsstellConnect.com to a USB memory. Insert the USB memory into the USB port of the Osstell IDx while it is on. The update will start automatically. Follow the instructions to finish.

### 21) Backup

A back-up of the Osstell IDx data is automatically saved if you are connected to WiFi or a fixed network. It can also be saved to a USB Memory stick, see section 8, Preparation and Set-up, Settings section, Net service and Net backup. Choose "export" to create a backup. Choose "import" to restore an existing backup.

### 22) Technical information

#### **Technical description**

Osstell IDx is CE-marked according to MDD in Europe (Class IIa, internally powered, type BF applied parts. Not AP or APG equipment, not protected against ingress of water).

Osstell IDx is in accordance with applicable parts of IEC 60601-1/ANSI/AAMI ES 60601-1.

The symbols used, follow the European standard EN 60601-1 and ISO 15223 as far as possible.

#### Notes on electromagnetic compatibility (EMC)

Medical electrical equipment is subjected to particular precautions with regards to EMC and must be put into operation in accordance with the EMC notes included below:

Osstell guarantees the compliance of the device with the EMC requirements only when used with original accessories and spare parts. The use of other accessories / other spare parts can lead to an increased emission of electromagnetic interference or to a reduced resistance against electromagnetic interference.

#### Battery and charger

The instrument contains a lithium-ion battery.

The instrument can only be charged with the Osstell IDx power supply. Battery charging is indicated by a blue LED surrounding the ON/OFF button. The battery symbol displays the status of the battery. A filled symbol indicates batteries are at 100% capacity.

The power supply operates with mains voltage from 100 to 240 VAC, 50-60 Hz. Thus it can be used in most countries in the world using the appropriate plug adapter. Use only in dry, indoor locations.

To disconnect the equipment from mains supply, disconnect the charger from the main socket outlet.

#### Accuracy

The Osstell IDx instrument has an ISQ accuracy/resolution of +/- 1 ISQ. When the SmartPeg is attached to an implant, the ISQ value can vary up to 2 ISQ depending on SmartPeg attachment torque.

#### Power, weight & size:

Rated power:	12 VDC
Charger type:	FW8001M/12
Instrument size:	203 x 163 x 72 mm
Package size:	270 x 243 x 105 mm
Instrument weight:	0.85 kg
Gross weight:	2.1 kg

#### Environmental conditions during transport:

Temperature:	-40 °C to +70 °C	
Relative humidity:	10% to 95%	
Pressure:	500 hPa to 1060 hPa	

#### Environmental conditions during use:

Temperature:	+10 °C to +35 °C	
Relative humidity:	30% to 75%	
Pressure:	700 hPa to 1060 hPa	
IP class:	IP20	

# 23) Troubleshooting

Problem	Cause	Possible Solution
No measurement or unexpected value	Re-Used SmartPeg.	The SmartPegs are disposable and should only be used for one or multiple measurements at one treatment session, for a single patient only. Repeated re-use may result in false readings due to wear and tear of the soft aluminium SmartPeg threads.
	Wrong SmartPeg type selected for the Implant.	See SmartPeg reference list, osstell.com/smartpegguide.
	Bone or soft tissue in between SmartPeg and Implant.	Make sure to clean the Implant prosthetic connection before attaching the SmartPeg.
	Electromagnetic interference. <b>(Noise)</b>	Remove the source of electromagnetic interference.
	Instrument tip is held too far away from the SmartPeg.	Normally it is sufficient to hold the instrument tip 3-5 mm away from the SmartPeg, but in some cases as close as 1 mm is necessary.
Unit is not charging when charger is connected	Wrong charger used.	Only use Osstell IDx with charger type FW8001M/12.
Charger does not fit the electrical outlets	Wrong Osstell plug adapter.	Use correct Osstell plug adapter.
Instrument does not start	Uncharged battery.	Charge the Osstell IDx.
Difficult to measure in an exact recommended direction	No space, e.g. due to adjacent teeth.	Try to measure at a slightly different angle.
Difficulties attaching the SmartPeg	Wrong SmartPeg.	Ensure that the SmartPeg is compatible with the implant system. See <b>osstell.com/smartpegguide</b> .
The instrument asks for registration	The Osstell IDx has not been registered.	Register the instrument at OsstellConnect.com Follow the step by step instructions on <b>osstell.com/get-started.</b>
The tooth view is not shown	Classic mode is activated.	Deactivate classic mode, see SETTINGS section.
# 24) Service and Support

In the event of an instrument malfunction contact the local sales representative or distributor for further instructions.

# 25) Waste and Disposal

The Osstell IDx (instrument, probe and charger) should be recycled as electrical equipment. SmartPegs should be recycled as metal. Whenever possible, the battery should be disposed in a discharged state to avoid heat generation through inadvertent short-circuiting.

Follow your local and country-specific laws, directives, standards and guidelines for disposal.

- Waste electrical equipment
- Accessories and spare parts
- Packaging

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