

SmartPeg Reference Guide

Valid from June 2016 | 25016-61

Implant Company	Model	Size	SmartPeg	Art. No
A.B. DENTAL DEVICES	Screw Type, Conical Groovy	3.75	Type 32	100440
	Screw Type, Conical Groovy	4.5	Type 27	100431
ADDENTA	Bifasico	4.0	Type 14	100374
ADIN DENTAL IMPLANT SYSTEM	One	3.3	Type 17	100390
	One	4.2	Type 17	100390
	Swell	3.3	Type 49	100480
	Swell	3.75	Type 49	100480
	Swell	4.2	Type 49	100480
	Swell	5.0	Type 49	100480
	Swell	6.0	Type 49	100480
	Touareg CloseFit	NP	Type 62	100545
	Touareg CloseFit	RP	Type 60	100543
	Touareg CloseFit	WP	Type 49	100480
	Touareg-S	4.2	Type 49	100480
Touareg-S	5.0	Type 49	100480	
AGS MEDIKAL	Implance Dental Implant System	Tissue Level	Type 06	100378
	Implance Dental Implant System	Bone Level	Type 06	100378
	Implance Dental Implant System	Mini Type	Type 38	100455
ALPHA BIO IMPLANT	ATI		Type 26	100425
	ATID	3.75	Type 26	100425
	NICE 3.2 Conical connection	3.2	Type 63	100554
	SPI		Type 26	100425
ALPHA-GATE	Bioactive	3.75	Type 35	100446
ALPHATECH		3.4	Type 46	100477
		3.8	Type 46	100477
		4.3	Type 46	100477
		5.0	Type 46	100477
		6.0	Type 46	100477
		7.0	Type 46	100477
ANKYLOS	C/X	A	Type 16	100388
	C/X	B	Type 16	100388
	C/X	C	Type 16	100388
ANTHOGYR	Anthofit External Hexagon HE		Type 01	100353
	Anthofit Internal Octagon OI		Type 18	100398
	Axiom PX	3.4	Type 38	100455

Implant Company	Model	Size	SmartPeg	Art. No
ANTHOGYR (Continued)	Axiom PX	4.0	Type 38	100455
	Axiom PX	4.6	Type 38	100455
	Axiom PX	5.2	Type 38	100455
	Axiom REG	3.4	Type 38	100455
	Axiom REG	4.0	Type 38	100455
	Axiom REG	4.6	Type 38	100455
	Axiom REG	5.2	Type 38	100455
	Ossfit		Type 04	100350
ARDS IMPLANTS	Classic	All	Type 27	100431
	Premium	All	Type 27	100431
	Smart	All	Type 27	100431
ASTRA TECH (Please see DENTSPLY)				
AVINENT IMPLANT SYSTEM	Coral External	3.5	Type 30	100436
	Coral External	4.1	Type 01	100353
	Coral External	5.1	Type 03	100357
	Ocean External	3.5	Type 30	100436
	Ocean External	4.1	Type 01	100353
	Internal Hexagon Coral	3.5	Type 32	100440
	Internal Hexagon Coral	4.1	Type 32	100440
	Internal Hexagon Ocean	3.5	Type 32	100440
	Internal Hexagon Ocean	4.1	Type 32	100440
AXIS	Ceramic	4.0	Type 38	100455
BEGO IMPLANT SYSTEMS	Semados	RI 3.75	Type 26	100425
	Semados	RI 4.1	Type 26	100425
	Semados	RI 4.5	Type 26	100425
	Semados	RI 5.5	Type 26	100425
	Semados	S 3.25	Type 26	100425
	Semados	S 3.75	Type 26	100425
	Semados	S 4.1	Type 26	100425
	Semados	S 4.5	Type 26	100425
	Semados	S 5.5	Type 26	100425
	Semados Mini	2.7	Type 02	100355
	Semados Mini	2.9	Type 02	100355
	Semados Mini	3.1	Type 02	100355
	Semados RS/R SX	3.0	Type 50	100481
	Semados RS/R SX	3.75	Type 26	100425
	Semados RS/R SX	4.1	Type 26	100425
	Semados RS/R SX	4.5	Type 26	100425
	Semados RS/R SX	5.5	Type 26	100425
BIODENTA	Tissue level	NP	Type 21	100402
	Tissue level	RP	Type 47	100478

Implant Company	Model	Size	SmartPeg	Art. No
BIODENTA (Continued)	Tissue level	WP	Type 47	100478
BIOHORIZONS	External	3.5	Type 30	100436
	External	4.0	Type 01	100353
	External	5.0	Type 01	100353
	External	6.0	Type 01	100353
	Internal	Blue platform	Type 33	100442
	Internal	Green platform	Type 27	100431
	Internal	Yellow platform	Type 32	100440
	Internal Laser Lok	3.0	Type 26	100425
	Single Stage	Blue platform	Type 33	100442
	Single Stage	Green platform	Type 27	100431
	Single Stage	Yellow platform	Type 32	100440
	Tapered Internal	Blue platform	Type 33	100442
	Tapered Internal	Green platform	Type 27	100431
	Tapered Internal	Yellow platform	Type 32	100440
BIOMET 3i	T3 Platform Switched	Purple platform	Type 45	100472
	T3 Platform Switched	Blue platform	Type 15	100386
	T3 Platform Switched	Yellow platform	Type 15	100386
	T3 Non-Platform Switched	Purple platform	Type 45	100472
	T3 Non-Platform Switched	Blue platform	Type 15	100386
	T3 Non-Platform Switched	Yellow platform	Type 15	100386
	T3 Non-Platform Switched	Green platform	Type 15	100386
	NanoTite External Hex	Purple platform	Type 01	100353
	NanoTite External Hex	Blue platform	Type 01	100353
	NanoTite External Hex	Yellow platform	Type 01	100353
	NanoTite External Hex	Green platform	Type 01	100353
	NanoTite External Hex	XP 3/4	Type 01	100353
	NanoTite External Hex	XP 4/5	Type 01	100353
	NanoTite External Hex	XP 5/6	Type 01	100353
	NanoTite External Hex MicroMiniplant	Purple platform	Type 01	100353
	NanoTite Tapered Certain Prevail	Purple platform	Type 45	100472
	NanoTite Tapered Certain Prevail	Blue platform	Type 15	100386
	NanoTite Tapered Certain Prevail	Yellow platform	Type 15	100386
	NanoTite Tapered Certain	Purple platform	Type 45	100472
	NanoTite Tapered Certain	Blue platform	Type 15	100386
	NanoTite Tapered Certain	Yellow platform	Type 15	100386
	NanoTite Tapered Certain	Green platform	Type 15	100386
	NanoTite Certain Microminiplant	3.25	Type 45	100472
	OSSEOTITE External Hex	Purple platform	Type 01	100353
	OSSEOTITE External Hex	Blue platform	Type 01	100353
	OSSEOTITE External Hex	Yellow platform	Type 01	100353
	OSSEOTITE External Hex	Green platform	Type 01	100353
	OSSEOTITE External Hex	XP 3/4	Type 01	100353
	OSSEOTITE External Hex	XP 4/5	Type 01	100353
	OSSEOTITE External Hex	XP 5/6	Type 01	100353

Implant Company	Model	Size	SmartPeg	Art. No
BIOMET 3i (Continued)	OSSEOTITE External Hex MicroMiniplant	Purple platform	Type 01	100353
	OSSEOTITE Tapered Certain Prevail	Purple platform	Type 45	100472
	OSSEOTITE Tapered Certain Prevail	Blue platform	Type 15	100386
	OSSEOTITE Tapered Certain Prevail	Yellow platform	Type 15	100386
	OSSEOTITE Tapered Certain	Purple platform	Type 45	100472
	OSSEOTITE Tapered Certain	Blue platform	Type 15	100386
	OSSEOTITE Tapered Certain	Yellow platform	Type 15	100386
	OSSEOTITE Tapered Certain	Green platform	Type 15	100386
	OSSEOTITE Certain	XP 4/5	Type 15	100386
	OSSEOTITE Certain	XP 5/6	Type 15	100386
	OSSEOTITE Certain Microminiplant	3.25	Type 45	100472
BIOTEC SRL	EN	All	Type 30	100436
	ER	All	Type 01	100353
	IM	All	Type 27	100431
	IR	All	Type 32	100440
	KR	3.05 platf	Type 38	100455
	KW	3.65 platf	Type 32	100440
	SR	All	Type 06	100378
	TWO	All	Type 32	100440
BONE SYSTEM	2P	5.0	Type 36	100447
BRAINBASE		C	Type 47	100478
	Mytis Arrow Implants	B	Type 04	100350
	Mytis Arrow Implants	BW	Type 04	100350
	Mytis Arrow Implants	E	Type 04	100350
	Mytis Arrow Implants	EW	Type 04	100350
BREDED	blueSKY		Type 49	100480
BTI BIOTECHNOLOGY INSTITUTE	Externa	Universal	Type 01	100353
	Externa	Universal plus	Type 01	100353
	Interna	Universal	Type 27	100431
	Interna	Universal plus	Type 27	100431
	Interna	Wide	Type 27	100431
	Tiny	2.5	Type 09	100368
	Tiny	3.75	Type 09	100368
BTLOCK	CV1	3.75	Type 49	100480
	CV1	4.5	Type 49	100480
	CV1	5.5	Type 49	100480
	CV1	6.5	Type 49	100480
	CV2	3.75	Type 49	100480
	CV2	4.5	Type 49	100480
	CV2	5.5	Type 49	100480
	CV2	6.5	Type 49	100480

Implant Company	Model	Size	SmartPeg	Art. No
BTLOCK (Continued)	CV3	3.75	Type 49	100480
	CV3	4.5	Type 49	100480
	CV3	5.5	Type 49	100480
	CV3	6.5	Type 49	100480
CAMLOG	Camlog Cylinder-line	3.3	Type 28	100432
	Camlog Cylinder-line	3.8	Type 23	100411
	Camlog Cylinder-line	4.3	Type 23	100411
	Camlog Cylinder-line	5.0	Type 24	100412
	Camlog Cylinder-line	6.0	Type 24	100412
	Camlog Root-line Promote	3.8	Type 23	100411
	Camlog Root-line Promote	4.3	Type 23	100411
	Camlog Root-line Promote	5.0	Type 24	100412
	Camlog Root-line Promote	6.0	Type 24	100412
	Camlog Screw-line	3.3	Type 28	100432
	Camlog Screw-line	3.8	Type 23	100411
	Camlog Screw-line	4.3	Type 23	100411
	Camlog Screw-line	5.0	Type 24	100412
	Camlog Screw-line	6.0	Type 24	100412
	Conelog Screw-line	3.3	Type 65	100557
	Conelog Screw-line	3.8	Type 58	100541
	Conelog Screw-line	4.3	Type 58	100541
	Conelog Screw-line	5.0	Type 07	100380
	iSy (direct to implant)	All	Type 57	100497
	iSy (to implant base)	All	Type 35	100446
CHANG GUNG	ComMed	3.5	Type 50	100481
	ComMed	4.0	Type 50	100481
	ComMed	4.5	Type 07	100380
	ComMed	5.0	Type 07	100380
	ComMed	5.5	Type 07	100380
	ComMed	6.0	Type 07	100380
CHAORUM	External	3.3	Type 01	100353
	Internal	4.0	Type 47	100478
	Internal	4.8	Type 47	100478
COCHLEAR	Baha BI 300	3.0	Type 30	100436
	Baha BI 300	4.0	Type 30	100436
	Baha BIA 300	6.0	Type 55	100487
	Baha BIA 300	9.0	Type 55	100487
	Baha BIA 400	6.0	Type 55	100487
	Baha BIA 400	8.0	Type 55	100487
	Baha BIA 400	10.0	Type 55	100487
	Baha BIA 400	12.0	Type 55	100487
	Baha BIA 400	14.0	Type 55	100487
	Baha & Vistafix Flange Fixture Hexagonal Connection	3.0	Type 09	100368

Implant Company	Model	Size	SmartPeg	Art. No
COCHLEAR (Continued)	Baha & Vistafix Flange Fixture Hexagonal Connection	4.0	Type 09	100368
	Vistafix VXI 300	3.0	Type 30	100436
	Vistafix VXI 300	4.0	Type 30	100436
COMPRESS		4.0	Type 02	100355
CORTEX	Classix	3.3	Type 26	100425
	Classix	3.8	Type 26	100425
	Classix	4.2	Type 26	100425
	Classix	5.0	Type 26	100425
	Classix	6.0	Type 26	100425
	Dynamix	3.3	Type 26	100425
	Dynamix	3.8	Type 26	100425
	Dynamix	4.2	Type 26	100425
	Dynamix	5.0	Type 26	100425
	Dynamix	6.0	Type 26	100425
	Saturn	3.3	Type 26	100425
	Saturn	3.8	Type 26	100425
	Saturn	4.2	Type 26	100425
	Saturn	5.0	Type 26	100425
Saturn	6.0	Type 26	100425	
COWELLMEDI	Atlas External	3.3	Type 02	100355
	Atlas External	4.1	Type 01	100353
	Atlas External	5.0	Type 01	100353
	Atlas Submerged	4.0	Type 21	100402
C-TECH	EL	3.5	Type 26	100425
	EL	4.3	Type 26	100425
	EL	5.1	Type 26	100425
DENTAL TECH		WP	Type 47	100478
DENTAURUM	Tiologic	S	Type 45	100472
	Tiologic	M	Type 45	100472
	Tiologic	L	Type 45	100472
DENTEGRIS	Straight & Tapered	3.3	Type 26	100425
	Straight & Tapered	3.75	Type 32	100440
	Straight & Tapered	4.5	Type 32	100440
	Straight & Tapered	4.75	Type 32	100440
	Straight & Tapered	5.5	Type 47	100478
DENTIS DENTAL	Internal	3.7	Type 04	100350
	Internal	4.3	Type 04	100350
	Internal	4.8	Type 04	100350
	OneQ S-clean Narrow (SLA surface/Submerged)	3.0	Type 22	100404

Implant Company		Model	Size	SmartPeg	Art. No
DENTIS DENTAL (Continued)		OneQ S-clean Narrow (SLA surface/Submerged)	3.3	Type 22	100404
		OneQ S-clean Regular (SLA surface/Submerged)	3.7	Type 11	100372
		OneQ S-clean Regular (SLA surface/Submerged)	4.2	Type 11	100372
		OneQ S-clean Regular (SLA surface/Submerged)	4.7	Type 11	100372
		OneQ S-clean Regular (SLA surface/Submerged)	5.2	Type 11	100372
		OneQ S-clean Wide (SLA surface/Submerged)	6.0	Type 11	100372
		OneQ S-clean Wide (SLA surface/Submerged)	7.0	Type 11	100372
		OneQ S-clean Wide (SLA surface/Submerged)	8.0	Type 11	100372
		Save Submerged	5.5	Type 07	100380
		Save Submerged	6.0	Type 07	100380
		S-clean Mini (Submerged)	3.7	Type 11	100372
		S-clean Regular (Submerged)	4.1	Type 11	100372
		S-clean Regular (Submerged)	4.3	Type 11	100372
		S-clean Wide (Submerged)	4.8	Type 11	100372
	DENTIUM		Genoss	3.0	Type 63
		Genoss	3.6	Type 21	100402
		Genoss	4.3	Type 21	100402
		Genoss	5.0	Type 21	100402
		Implantium	3.4 (3.6)	Type 06	100378
		Implantium	3.8 (4.0)	Type 06	100378
		Implantium	4.3 (4.5)	Type 06	100378
		Simple Line	3.8	Type 07	100380
		Simple Line	4.3	Type 07	100380
		Simple Line	4.8	Type 07	100380
	Super Line	4.0	Type 07	100380	
	Super Line	7.0	Type 07	100380	
DENTSPLY IMPLANTS		ASTRA TECH Implant system MicroThread	3.5	Type 05	100376
		ASTRA TECH Implant system MicroThread	4.0	Type 06	100378
		ASTRA TECH Implant system MicroThread	4.5	Type 07	100380
		ASTRA TECH Implant system MicroThread	5.0	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed Yellow	3.0S	Type 50	100481
		ASTRA TECH Implant system OsseoSpeed Aqua	3.5S	Type 38	100455
		ASTRA TECH Implant system OsseoSpeed Aqua	4.0S	Type 38	100455
		ASTRA TECH Implant system OsseoSpeed Lilac	4.5	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed Lilac	5.0	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed Lilac	5.0S	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed TX Yellow	3.0S	Type 50	100481
		ASTRA TECH Implant system OsseoSpeed TX Aqua	3.5S	Type 38	100455
		ASTRA TECH Implant system OsseoSpeed TX Aqua	4.0S	Type 38	100455
		ASTRA TECH Implant system OsseoSpeed TX Lilac	4.5	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed TX Lilac	5.0	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed TX Lilac	5.0S	Type 07	100380
		ASTRA TECH Implant system OsseoSpeed TX Profile	4.5	Type 06	100378
		ASTRA TECH Implant system OsseoSpeed TX Profile	5.0	Type 06	100378
		ASTRA TECH Implant system OsseoSpeed TX Profile	5.0S	Type 06	100378

Implant Company	Model	Size	SmartPeg	Art. No
DENTSPLY IMPLANTS (Contin.)	ASTRA TECH Implant system OsseoSpeed EV	4.2C Yellow	Type 49	100480
	ASTRA TECH Implant system OsseoSpeed EV	4.8C Yellow	Type 47	100478
	ASTRA TECH Implant system OsseoSpeed EV	3.0S Green	Type 50	100481
	ASTRA TECH Implant system OsseoSpeed EV	3.6S Purple	Type 38	100455
	ASTRA TECH Implant system OsseoSpeed EV	4.2S Yellow	Type 49	100480
	ASTRA TECH Implant system OsseoSpeed EV	4.8S Blue	Type 47	100478
	ASTRA TECH Implant system OsseoSpeed EV	5.4S Brown	Type 40	100462
	ASTRA TECH Implant system OsseoSpeed Profile EV	4.2 PC	Type 10	100370
	ASTRA TECH Implant system OsseoSpeed Profile EV	4.8 PC	Type 32	100440
	ASTRA TECH Implant system OsseoSpeed Profile EV	4.2 PS	Type 10	100370
	ASTRA TECH Implant system OsseoSpeed Profile EV	4.8 PS	Type 32	100440
	ANKYLOS C/X Implants	A	Type 16	100388
	ANKYLOS C/X Implants	B	Type 16	100388
	ANKYLOS C/X Implants	C	Type 16	100388
	XiVe	3.0	Type 22	100404
	XiVe	3.4	Type 21	100402
	XiVe	3.8	Type 45	100472
	XiVe	4.5	Type 45	100472
	XiVe	5.5	Type 45	100472
	XiVe	TG	Type 08	100366
DIO	External FSN	3.3	Type 02	100355
	External FSN	3.75	Type 01	100353
	External FTN	5.0	Type 03	100357
	SM Submerged	3.8	Type 49	100480
	SM Submerged	4.5	Type 47	100478
	SM Submerged	5.3	Type 47	100478
DTI		3.5	Type 57	100497
		4.2	Type 07	100380
DYNA DENTAL	Dyna Helix DC	3.2	Type 27	100431
	Dyna Helix DC	3.6	Type 27	100431
	Dyna Helix DC	4.2	Type 27	100431
	Dyna Helix DC	5.0	Type 27	100431
	Dyna Helix ST	3.6	Type 27	100431
	Dyna Helix ST	4.2	Type 27	100431
	Dyna Helix ST	5.0	Type 27	100431
ECKERMANN	All Spiral		Type 17	100390
	Hexagono Ext Hex		Type 01	100353
	Triplo		Type 49	100480
EMFILS	Colosso IH ECI Evolution	3.3 x 16	Type 11	100372
	Colosso IH ECI Evolution	4.0 x 7	Type 11	100372
	Colosso IH ECI Evolution	5.0 x 16	Type 11	100372
	ICM-4010	4.0 x 10	Type 53	100484

Implant Company	Model	Size	SmartPeg	Art. No
EMFILS (Continued)	PAT5-4010	4.0 x 10	Type 53	100484
ENTIFIC	Flange Fixture		Type 01	100353
EUROTEKNIKA	Aesthetica+2	3.6/4.2	Type 61	100544
	Aesthetica+2	3.6/4.8	Type 61	100544
	Aesthetica+2	4.1/4.2	Type 61	100544
	Aesthetica+2	4.1/4.8	Type 61	100544
	Aesthetica+2	4.8/4.8	Type 61	100544
	Aesthetica+2	4.8/6.5	Type 61	100544
	Naturall	3.5	Type 38	100455
	Naturall	4.0	Type 06	100378
	Naturall	4.5	Type 06	100378
	Naturall	5.0	Type 06	100378
	Uneva+	3.6	Type 01	100353
	Uneva+	4.1	Type 01	100353
	Uneva+	4.8	Type 01	100353
	Uneva+	6.0	Type 01	100353
FAIR IMPLANT	FairTwo	3.5	Type 53	100484
	FairTwo	4.2	Type 53	100484
	FairTwo	5.0	Type 06	100378
	FairTwo	6.0	Type 06	100378
GC IMPLANT	Aadva	Narrow	Type 38	100455
	Aadva	Regular	Type 58	100541
	Aadva	Wide	Type 58	100541
	External Plus	3.0	Type 21	100402
	External Plus	3.4	Type 21	100402
	External Plus	3.8	Type 01	100353
	External Plus	5.0	Type 01	100353
	GENESIO		Type 21	100402
	Internal Plus	3.0	Type 21	100402
	Internal Plus	3.4	Type 21	100402
	SETIO		Type 01	100353
GHIMAS	PFON EVO SERIES	3.6	Type 38	100455
	PFON EVO SERIES	4.0	Type 38	100455
	PFON EVO SERIES	5.0	Type 07	100380
HIOSSEN	ETIII	4.0	Type 47	100478
	ETIII	4.5	Type 47	100478
	ETIII	5.0	Type 47	100478
	ETIII	Mini 3.5	Type 05	100376
	ETIII	UW 6.0	Type 06	100378
	ETIII	UW 7.0	Type 06	100378
	HGII	Mini	Type 05	100376

Implant Company	Model	Size	SmartPeg	Art. No
HIOSSSEN (Continued)	HGII	Regular	Type 06	100378
	HGII	Wide	Type 06	100378
		Ultra Wide	Type 06	100378
	HGIII	Mini	Type 05	100376
	HSII	Regular	Type 04	100350
	HSII	Wide	Type 04	100350
IBS INNO BIO SURG	Magicore	5.5	Type 38	100455
IDI EVOLUTION	Vite	3.75	Type 18	100398
	Vite	4.0	Type 18	100398
IMMEDIATELOAD	IM40 Immediate Conico	4.0	Type 49	100480
IMPLANT DIRECT	Legacy3	3.2	Type 38	100455
	ScrewPlant	3.7	Type 26	100425
	ScrewPlant	4.7	Type 27	100431
	ScrewPlant	5.7	Type 27	100431
IMTEC	Hexed-head	3.75	Type 01	100353
	Hexed-head	4.75	Type 01	100353
INTRA-LOCK SYSTEM	Internal	SQ	Type 06	100378
	Internal	Standard	Type 47	100478
	Internal	Wide	Type 47	100478
ITERUM	Conical	3.75	Type 26	100425
JDENTALCARE	JDEvolution	3.2	Type 38	100455
	JDEvolution	3.7	Type 32	100440
	JDEvolution	4.3	Type 32	100440
	JDEvolution	5.0	Type 32	100440
JMM	POI 2 piece	3.2	Type 18	100398
	POI 2 piece	3.7	Type 18	100398
	POI 2 piece	4.2	Type 18	100398
	POIEX	3.7	Type 15	100386
	POIEX	4.2	Type 49	100480
	POIEX	4.7	Type 47	100478
	POIEX	5.2	Type 47	100478
KENTEC	SB1	3.3	Type 06	100378
	SB2	4.2	Type 06	100378
KEYSTONE DENTAL	Genesis	3.8	Type 27	100431
	Genesis	4.5	Type 33	100442
	Genesis	5.5	Type 33	100442

Implant Company	Model	Size	SmartPeg	Art. No
KEYSTONE DENTAL (Continued)	Genesis	6.5	Type 33	100442
	Prima Connex	SD 3.5	Type 27	100431
	Prima Connex	RD 4.1	Type 33	100442
	Prima Connex	WD 5.0	Type 33	100442
	Renova SDI	3.75	Type 26	100425
	Renova RDI	4.5	Type 27	100431
	Restore	3.3	Type 30	100436
	Restore	3.75	Type 01	100353
	Restore	4.0	Type 01	100353
	Restore	5.0	Type 03	100357
	Restore	6.0	Type 03	100357
KLOCKNER IMPLANT SYSTEM	MV VEGA Implant	3.0	Type 57	100497
	NV VEGA Implant	3.5	Type 57	100497
	RV VEGA Implant	4.0	Type 26	100425
	RV VEGA Implant	4.5	Type 26	100425
	Essential Cone	3.5	Type 49	100480
	Essential Cone	4.0	Type 49	100480
	Essential Cone	4.5	Type 49	100480
	Essential Cone	4.8	Type 49	100480
	KL Implant NP	3.3	Type 02	100355
	KL Implant NP	3.7	Type 01	100353
	KL Implant NP	4.2	Type 01	100353
	KL Implant WP	4.7	Type 03	100357
	S3M		Type 30	100436
	SK2		Type 30	100436
KYOCERA	FINESIA BL	NP	Type 62	100545
	FINESIA BL	RP	Type 27	100431
	FINESIA BL	WP	Type 04	100350
	FINESIA TL	RP	Type 04	100350
	FINESIA TL	WP	Type 04	100350
LASAK	BioniQ Q-lock	Regular	Type 57	100497
	Impladent	2.9	Type 02	100355
	Impladent	3.7	Type 36	100447
	Impladent	5.0	Type 37	100448
LEADER	Implus EH	3.3	Type 01	100353
	Implus EH	3.75	Type 01	100353
	Implus EH	5.0	Type 01	100353
	Implus IH	3.3	Type 26	100425
	Tixos IH	3.75	Type 26	100425
	Tixos IH	4.5	Type 26	100425
	Tixos IH	5.5	Type 26	100425
	Tixos ITMC	3.75	Type 57	100497
	Tixos ITMC	4.5	Type 57	100497

Implant Company	Model	Size	SmartPeg	Art. No
LEADER (Continued)	Tixos ITMC	5.5	Type 57	100497
MEDENTIKA	M-Implant®	3.5	Type 38	100455
	M-Implant®	4.0	Type 38	100455
	M-Implant®	4.5	Type 38	100455
	M-Implant®	5.0	Type 38	100455
MEDICAL INSTINCT	Bone Trust	3.4	Type 44	100466
	Bone Trust	4.0	Type 44	100466
	Bone Trust	5.0	Type 44	100466
	Bone Trust Plus Cone	3.4	Type 44	100466
	Bone Trust Plus Cone	4.0	Type 44	100466
MEDIGMA TECHNOLOGIES	Fix 1		Type 17	100390
	Fix-A-Dent		Type 49	100480
	FixTite		Type 49	100480
	FixTite-S		Type 49	100480
MEGAGEN	AnyOne	3.5	Type 07	100380
	AnyOne	4.0	Type 07	100380
	AnyOne External	4.5	Type 04	100350
	AnyRidge		Type 27	100431
	EZ Plus Internal	Small	Type 38	100455
	EZ Plus Internal	Regular	Type 07	100380
	EZ Plus Internal	Wide	Type 07	100380
	EZ Plus Internal	4.0	Type 07	100380
	EZ Plus Internal	4.5	Type 07	100380
	EZ Plus Internal	5.0	Type 07	100380
	Mini Internal	3.0	Type 50	100481
	Resque External	7.0	Type 03	100357
MICRODENT	Duraplant 2.2	4.1	Type 10	100370
	Microdent system		Type 01	100353
	MK	5.6	Type 03	100357
	MS	3.5	Type 30	100436
	MST	3.5	Type 30	100436
	Trylogic	3.8	Type 27	100431
	Trylogic	4.3	Type 27	100431
	Universal	4.1	Type 01	100353
	Universal	5.1	Type 01	100353
MIS	C1 (conical connection)	Standard 3.75	Type 49	100480
	Internal Hexagon (seven, biocom)	3.75	Type 32	100440
	Internal Octagon (mistral)	4.75	Type 04	100350
	Seven Narrow	3.3	Type 63	100554

Implant Company	Model	Size	SmartPeg	Art. No
MOZO-GRAU	MG-Inhex	Std	Type 38	100455
	MG-Inhex	Maxi	Type 07	100380
	MG-Osseous	Mini	Type 01	100353
	MG-Osseous	Std	Type 01	100353
	MG-Osseous	Maxi	Type 03	100357
NEMRIS	Aesthura Classic	3.75	Type 35	100446
	Aesthura Classic	4.75	Type 35	100446
NEOBIOTECH	IS Fixture (RBM Surface)	4.0	Type 05	100376
	IS Fixture (RBM Surface)	4.5	Type 05	100376
	IS Fixture (RBM Surface)	4.8	Type 05	100376
	IS S-Wide Fixture (RBM Surface)	5.5	Type 05	100376
	IS S-Wide Fixture (RBM Surface)	6.0	Type 05	100376
	IS S-Wide Fixture (RBM Surface)	7.0	Type 05	100376
	IS-II Fixture (RBM Surface)	3.5	Type 05	100376
	IS-II Fixture (RBM Surface)	4.0	Type 05	100376
	IS-II Fixture (RBM Surface)	4.5	Type 05	100376
	IS-II Fixture (RBM Surface)	5.0	Type 05	100376
	IS-II S-Wide Fixture (RBM Surface)	5.5	Type 05	100376
	IS-II S-Wide Fixture (RBM Surface)	6.0	Type 05	100376
	IS-II S-Wide Fixture (RBM Surface)	7.0	Type 05	100376
	IS-II active Fixture (S.L.A. Surface)	3.5	Type 05	100376
	IS-II active Fixture (S.L.A. Surface)	4.0	Type 05	100376
	IS-II active Fixture (S.L.A. Surface)	4.5	Type 05	100376
	IS-II active Fixture (S.L.A. Surface)	5.0	Type 05	100376
	IS-II active S-Wide Fixture (S.L.A. Surface)	5.5	Type 05	100376
	IS-II active S-Wide Fixture (S.L.A. Surface)	6.0	Type 05	100376
	IS-II active S-Wide Fixture (S.L.A. Surface)	7.0	Type 05	100376
	IS-III active Fixture	3.5	Type 05	100376
	IS-III active Fixture	4.0	Type 05	100376
	IS-III active Fixture	4.5	Type 05	100376
	IS-III active Fixture	5.0	Type 05	100376
	IS-III active S-Wide Fixture	5.5	Type 05	100376
	IS-III active S-Wide Fixture	6.0	Type 05	100376
	IS-III active S-Wide Fixture	7.0	Type 05	100376
	Neoplant External	4.0	Type 01	100353
	Neoplant External	5.0	Type 01	100353
	NEODENT	Alvim CM	3.5	Type 16
Alvim CM		4.3	Type 16	100388
Alvim CM		5.0	Type 16	100388
Alvim II Plus		3.75	Type 21	100402
Alvim II Plus		4.3	Type 21	100402
Alvim II Plus		5.0	Type 21	100402
Alvim Ti		4.3	Type 01	100353
Drive CM	4.3	Type 48	100479	

Implant Company	Model	Size	SmartPeg	Art. No
NEODENT (Continued)	Drive CM	5.0	Type 48	100479
	Titamax CM	3.5	Type 16	100388
	Titamax CM	3.75	Type 16	100388
	Titamax CM	4.0	Type 16	100388
	Titamax CM	5.0	Type 16	100388
	Titamax GT	3.7	Type 09	100368
	Titamax GT	4.0	Type 09	100368
	Titamax GT	4.5	Type 09	100368
	Titamax II Plus	3.75	Type 21	100402
	Titamax II Plus	5.0	Type 21	100402
	Titamax Ti	3.75	Type 01	100353
	Titamax Ti	4.0	Type 01	100353
	Titamax Ti	4.5	Type 03	100357
	Titamax Ti	5.0	Type 03	100357
NEOSS	Neoss Implant	3.25	Type 38	100455
	ProActive Straight	3.5	Type 18	100398
	ProActive Straight	4.0	Type 18	100398
	ProActive Straight	4.5	Type 18	100398
	ProActive Straight	5.0	Type 18	100398
	ProActive Straight	5.5	Type 18	100398
	ProActive Straight	6.0	Type 18	100398
	ProActive Tapered	3.5	Type 18	100398
	ProActive Tapered	4.0	Type 18	100398
	ProActive Tapered	4.5	Type 18	100398
	ProActive Tapered	5.0	Type 18	100398
	ProActive Tapered	5.5	Type 18	100398
	ProActive Tapered	6.0	Type 18	100398
NIPPON PISTON RING	1 Stage Screw RH	3.3	Type 18	100398
	1 Stage Screw RH	4.0	Type 18	100398
	1 Stage Screw RL	3.3	Type 18	100398
	1 Stage Screw RL	4.0	Type 18	100398
	1 Stage Screw WH	5.0	Type 18	100398
	1 Stage Screw WL	5.0	Type 18	100398
	2 Stage Cylinder	3.3	Type 18	100398
	2 Stage Cylinder	4.0	Type 18	100398
	2 Stage Cylinder	5.0	Type 18	100398
	2 Stage Screw	3.3	Type 18	100398
	2 Stage Screw	4.0	Type 18	100398
	2 Stage Screw	5.0	Type 18	100398
NOBEL BIOCARE	Brånemark system	NP	Type 02	100355
	Brånemark system	RP	Type 01	100353
	Brånemark system	WP	Type 03	100357
	Brånemark system Novum		Type 01	100353
	Brånemark system Shorty	RP	Type 01	100353

Implant Company	Model	Size	SmartPeg	Art. No
NOBEL BIOCARE (Continued)	Brånemark system Shorty	WP	Type 03	100357
	NobelActive	NP	Type 60	100543
	NobelActive	RP	Type 61	100544
	NobelActive	WP	Type 61	100544
	NobelActive	WP 5.5	Type 04	100350
	NobelActive	3.0	Type 62	100545
	NobelDirect	NP	Type 08	100366
	NobelDirect	RP	Type 09	100368
	NobelDirect	WP	Type 09	100368
	NobelParallel Conical Connection	NP	Type 60	100543
	NobelParallel Conical Connection	RP	Type 61	100544
	NobelParallel Conical Connection	WP	Type 61	100544
	NobelReplace Conical Connection	NP	Type 60	100543
	NobelReplace Conical Connection	RP	Type 61	100544
	NobelReplace Conical Connection PMC	NP	Type 60	100543
	NobelReplace Conical Connection PMC	RP	Type 61	100544
	NobelReplace Platform Shift	NP	Type 12	100359
	NobelReplace Platform Shift	RP	Type 13	100361
	NobelReplace Platform Shift	WP	Type 13	100361
	NobelReplace Straight Groovy	NP	Type 12	100359
	NobelReplace Straight Groovy	RP	Type 13	100361
	NobelReplace Straight Groovy	WP	Type 13	100361
	NobelReplace Tapered Groovy	NP	Type 12	100359
	NobelReplace Tapered Groovy	RP	Type 13	100361
	NobelReplace Tapered Groovy	WP	Type 13	100361
	NobelReplace Tapered Groovy	6.0	Type 40	100462
	NobelSpeedy Groovy	NP	Type 02	100355
	NobelSpeedy Groovy	RP	Type 01	100353
	NobelSpeedy Groovy	WP	Type 03	100357
	NobelSpeedy Replace	NP	Type 12	100359
	NobelSpeedy Replace	RP	Type 13	100361
	NobelSpeedy Replace	WP	Type 13	100361
	NobelSpeedy Replace	6.0	Type 40	100462
	NobelSpeedy Shorty	RP	Type 01	100353
	NobelSpeedy Shorty	WP	Type 03	100357
	Replace Select Straight	NP	Type 12	100359
	Replace Select Straight	RP	Type 13	100361
	Replace Select Straight	WP	Type 13	100361
	Replace Select Tapered	NP	Type 12	100359
	Replace Select Tapered	RP	Type 13	100361
	Replace Select Tapered	WP	Type 13	100361
	Replace Select Tapered	6.0	Type 40	100462
	Replace Select Tapered PMC	NP	Type 12	100359
	Replace Select Tapered PMC	RP	Type 13	100361
	Replace Select Tapered PMC	WP	Type 13	100361
	Replace Select Tapered PMC	6.0	Type 40	100462
	Replace Select TC	NP	Type 12	100359

Implant Company	Model	Size	SmartPeg	Art. No
NOBEL BIOCARE (Continued)	Replace Select TC	RP	Type 13	100361
NUCLEOSS	Nucleoss T3	3.4	Type 05	100376
	Nucleoss T3	3.8	Type 05	100376
	Nucleoss T3	4.2	Type 05	100376
	Nucleoss T3	5.0	Type 05	100376
	Nucleoss T4	3.4	Type 05	100376
	Nucleoss T4	3.8	Type 05	100376
	Nucleoss T4	4.2	Type 05	100376
	Nucleoss T5	5.0	Type 07	100380
OCO BIOMEDICAL	Enagage	3.2	Type 38	100455
	Enagage	3.7/3.5	Type 26	100425
	Enagage	4.7/4.5	Type 26	100425
	Enagage	5.7/4.5	Type 27	100431
	ERI	4.0	Type 07	100380
	ERI	5.0	Type 07	100380
	Macro	6.0	Type 47	100478
	TSI	3.25	Type 35	100446
	TSI	4.0	Type 07	100380
	TSI	5.0	Type 07	100380
OSSEOLINK		3.5	Type 21	100402
		4.0	Type 21	100402
		5.0	Type 21	100402
OSSTEM	GSII	Mini	Type 05	100376
	GSII	Regular	Type 06	100378
	GSII	Ultra Wide	Type 06	100378
	GSII	Wide	Type 06	100378
	GSIII	Mini	Type 05	100376
	GSIII	Regular	Type 06	100378
	SSI	Regular	Type 04	100350
	SSI	Wide	Type 04	100350
	SSII	Regular	Type 04	100350
	SSII	Wide	Type 04	100350
	SSIII	Regular	Type 04	100350
	SSIII	Wide	Type 04	100350
	TSIII CA	3.5	Type 05	100376
	TSIII SA	3.0	Type 11	100372
	TSIII SA	4.0	Type 47	100478
	TSIII SA	4.5	Type 47	100478
	TSIII SA	5.0	Type 47	100478
	TSIII SA Fixture – NoMount	4.5	Type 47	100478
	TSIII SA Ultra-Wide	6.0	Type 04	100350
	TSIII SA Ultra-Wide	7.0	Type 04	100350
	USII	Mini	Type 02	100355

Implant Company	Model	Size	SmartPeg	Art. No
OSSTEM (Continued)	USII	Regular	Type 01	100353
	USII	Ultra Wide	Type 03	100357
	USII	Wide	Type 03	100357
	USIII	Mini	Type 02	100355
	USIII	Regular	Type 01	100353
	USIII	Ultra Wide	Type 03	100357
	USIII	Wide	Type 03	100357
	USIV Conical	Mini	Type 08	100366
OSTEOPLANT	Two-piece		Type 10	100370
OSTEO-TI	Combi Dental Implant	3.25	Type 26	100425
	Combi Dental Implant	3.75	Type 26	100425
	Combi Dental Implant	4.5	Type 26	100425
	Combi Dental Implant	5.5	Type 26	100425
OTICON MEDICAL	Ponto Implant	3.0	Type 09	100368
	Ponto Implant	4.0	Type 09	100368
	Ponto Implant	Wide	Type 09	100368
PALTOP	PALTOP	NP	Type 38	100455
	PALTOP	SP	Type 49	100480
PHIBO	Aurea	NP	Type 63	100554
	Aurea	RP	Type 27	100431
	Aurea	WP	Type 64	100555
	Aurea EVO	NP	Type 63	100554
	Aurea EVO	RP	Type 26	100425
	Aurea EVO	WP	Type 26	100425
	BNT	S3	Type 01	100353
	BNT	S4	Type 01	100353
	BNT	S5	Type 01	100353
	TSA	S3	Type 38	100455
	TSA	S4	Type 27	100431
	TSA	S5	Type 27	100431
	TSA Advance	S3	Type 38	100455
	TSA Advance	S4	Type 27	100431
	TSA Advance	S5	Type 27	100431
	TSH	S3	Type 01	100353
	TSH	S4	Type 01	100353
	TSH	S5	Type 01	100353
P-I BRÄNEMARK PHILOSOPHY	Internal Grip Amplified	3.5	Type 63	100554
	Internal Grip Amplified	4.3	Type 51	100482
	Internal Grip Amplified	5.1	Type 51	100482
	Internal Grip External Hex	3.5	Type 38	100455
	Internal Grip External Hex	3.75	Type 18	100398

Implant Company	Model	Size	SmartPeg	Art. No
P-I BRÄNEMARK PHILOSOPHY	Internal Grip External Hex	4.1	Type 18	100398
(Continued)	Morse Taper	3.3	Type 21	100402
	Morse Taper	3.75	Type 21	100402
	Morse Taper	4.8	Type 21	100402
PLAN 1 HEALTH	GTB	3.3	Type 57	100497
	GTB	3.6	Type 57	100497
	GTB	4.3	Type 57	100497
	Hexcel EHC, External Hex		Type 01	100353
	Hexcel IHC, Internal Hex		Type 26	100425
	Hexcel IHT, Internal Hex		Type 26	100425
PLATON	EL Implant Narrow	3.7	Type 60	100543
	EL Implant Regular	5.1	Type 47	100478
	Type I		Type 04	100350
	Type II		Type 04	100350
	Type IV Pro		Type 47	100478
PRODENT ITALIA	Perform SM	3.6	Type 49	100480
	Perform SM	4.0	Type 49	100480
	Perform SM	4.5	Type 07	100380
	Perform SM	5.0	Type 07	100380
REVIOS		3.8	Type 26	100425
		4.3	Type 26	100425
		5.0	Type 26	100425
REX IMPLANTS	Rex TL	1.8	Type 04	100350
	Rex TL	2.9	Type 49	100480
SDS	SDS 2.0	3.8	Type 69	100589
SERVO-DENTAL	Trias	3.8	Type 06	100378
SHINHUNG	Luna	NP 3.5	Type 05	100376
	Luna	RP 4.0	Type 47	100478
	Luna	RP 4.5	Type 47	100478
	Luna	RP 5.0	Type 47	100478
	Sola	NP 3.5	Type 02	100355
	Sola	RP 4.0	Type 01	100353
	Sola	RP 4.5	Type 01	100353
	Sola	WP 5.0	Type 03	100357
	Stella	RP 4.0	Type 61	100544
	Stella	RP 4.5	Type 61	100544
	Stella	WP 4.5	Type 61	100544
	Stella	WP 5.0	Type 61	100544

Implant Company		Model	Size	SmartPeg	Art. No
SIC IMPLANT SYSTEM		Ace	3.4	Type 44	100466
		Ace	4.0	Type 44	100466
		Ace	5.0	Type 44	100466
S.I.N IMPLANT SYSTEM		External Hex	3.4	Type 01	100353
		External Hex	4.1	Type 01	100353
		External Hex	5.0	Type 01	100353
		Internal Hex	3.8	Type 15	100386
		Internal Hex	4.5	Type 15	100386
		Internal Hex	5.5	Type 15	100386
		Strong SW HI	4.5	Type 48	100479
		Strong SW Morse	4.5	Type 53	100484
SOUTHERN IMPLANTS		External Hex		Type 01	100353
		Extra-oral	3.75	Type 01	100353
		IT Connection		Type 04	100350
		MAX	7.0	Type 01	100353
		MAX	8.0	Type 01	100353
		MAX	9.0	Type 01	100353
		MAX-IT	7.0	Type 07	100380
		MAX-IT	8.0	Type 07	100380
		MAX-IT	9.0	Type 07	100380
		Tri-MAX	7.0	Type 24	100412
		Tri-MAX	8.0	Type 40	100462
		Tri-MAX	9.0	Type 40	100462
		Tri-Nex	3.5	Type 12	100359
		Tri-Nex	4.3	Type 13	100361
		Tri-Nex	5.0	Type 13	100361
STERNGOLD DENTAL		External Hexagon	RP 3.75	Type 01	100353
		External Hexagon	RP 4.0	Type 01	100353
		External Hexagon	RP 5.0	Type 01	100353
		Stern IC Internal Octagon	RN 4.1	Type 04	100350
STRAUMANN	Bone Level Implants Roxolid SLA Active NC		3.3	Type 53	100484
	Bone Level Implants Roxolid SLA Active RC		4.1	Type 54	100485
	Bone Level Implants Roxolid SLA Active RC		4.8	Type 54	100485
	Bone Level Implants Titanium SLA NC		3.3	Type 53	100484
	Bone Level Implants Titanium SLA RC		4.1	Type 54	100485
	Bone Level Implants Titanium SLA RC		4.8	Type 54	100485
	Bone Level Tapered Implants Roxolid SLA Active NC		3.3	Type 53	100484
	Bone Level Tapered Implants Roxolid SLA Active RC		4.1	Type 54	100485
	Bone Level Tapered Implants Roxolid SLA Active RC		4.8	Type 54	100485
	Extra Oral			Type 01	100353
	Narrow Neck			Type 17	100390
	Ortho			Type 09	100368
	Palatal			Type 09	100368

Implant Company	Model	Size	SmartPeg	Art. No
STRAUMANN (Continued)	Standard Implants Roxolid SLA Active RN	3.3	Type 04	100350
	Standard Implants Roxolid SLA Active RN	4.1	Type 04	100350
	Standard Implants Roxolid SLA Active RN	4.8	Type 04	100350
	Standard Implants Roxolid SLA Active WN	4.8	Type 04	100350
	Standard Implants Titanium SLA RN	3.3	Type 04	100350
	Standard Implants Titanium SLA RN	4.1	Type 04	100350
	Standard Implants Titanium SLA RN	4.8	Type 04	100350
	Standard Implants Titanium SLA WN	4.8	Type 04	100350
	Standard Plus Implants Roxolid SLA Active NNC	3.3	Type 68	100560
	Standard Plus Implants Roxolid SLA Active RN	3.3	Type 04	100350
	Standard Plus Implants Roxolid SLA Active RN	4.1	Type 04	100350
	Standard Plus Implants Roxolid SLA Active RN	4.8	Type 04	100350
	Standard Plus Implants Roxolid SLA Active WN	4.8	Type 04	100350
	Standard Plus Implants Titanium SLA RN	3.3	Type 04	100350
	Standard Plus Implants Titanium SLA RN	4.1	Type 04	100350
	Standard Plus Implants Titanium SLA RN	4.8	Type 04	100350
	Standard Plus Implants Titanium SLA WN	4.8	Type 04	100350
	Tapered Effect Implants Roxolid SL Active RN	3.3	Type 04	100350
	Tapered Effect Implants Roxolid SL Active RN	4.1	Type 04	100350
	Tapered Effect Implants Roxolid SL Active WN	4.8	Type 04	100350
	Tapered Effect Implants Titanium SLA RN	3.3	Type 04	100350
	Tapered Effect Implants Titanium SLA RN	4.1	Type 04	100350
	Tapered Effect Implants Titanium SLA WN	4.8	Type 04	100350
SWEDEN & MARTINA	EZ		Type 07	100380
	Global	3.8	Type 48	100479
	Global	4.3	Type 48	100479
	Global	4.8	Type 48	100479
	Global	5.5	Type 48	100479
	Kohno - new platform	3.8	Type 32	100440
	Kohno - old platform	3.8	Type 21	100402
	Kohno - new platform	4.25	Type 07	100380
	Kohno - old platform	4.25	Type 07	100380
	Kohno - new platform	5.0	Type 07	100380
	Kohno - old platform	5.0	Type 07	100380
	Premium - new platform	3.3	Type 32	100440
	Premium - old platform	3.3	Type 21	100402
	Premium - new platform	3.8	Type 32	100440
	Premium - old platform	3.8	Type 21	100402
	Premium - new platform	4.25	Type 07	100380
	Premium - old platform	4.25	Type 07	100380
	Premium - new platform	5.0	Type 07	100380
	Premium - old platform	5.0	Type 07	100380
SYBRON	Endopore External	3.5	Type 30	100436
	Endopore External	4.1	Type 01	100353
	Endopore External	5.0	Type 01	100353

Implant Company	Model	Size	SmartPeg	Art. No
SYBRON (Continued)	Entegra External	3.25	Type 01	100353
	Entegra External	4.0	Type 01	100353
	Pitt-Easy	4.0	Type 24	100412
	Pitt-Easy	4.9	Type 24	100412
	SybronPRO XRT	3.3	Type 06	100378
	SybronPRO XRT	3.5	Type 06	100378
	SybronPRO XRT	4.1	Type 06	100378
TA-DENT		3.3	Type 38	100455
		4.8	Type 06	100378
		5.6	Type 49	100480
TBR	External Hexagon		Type 01	100353
	Internal Octagon		Type 18	100398
	M-implant	3.2	Type 21	100402
	M-implant	3.9	Type 21	100402
	M-implant	4.7	Type 21	100402
	Z1-M implant	3.2	Type 21	100402
	Z1-M implant	3.9	Type 21	100402
	Z1-M implant	4.7	Type 21	100402
TEKKA	Inkone	3.5	Type 49	100480
	Inkone	4.5	Type 49	100480
	Krestal	3.5	Type 06	100378
	Krestal	4.0	Type 06	100378
	Krestal	4.2	Type 06	100378
	Krestal	5.1	Type 06	100378
	Progress	3.4	Type 06	100378
	Progress	4.0	Type 06	100378
	T-quest	4.0	Type 49	100480
	T-quest	4.5	Type 49	100480
	Twinkone	3.5	Type 17	100390
TFI SYSTEM	TFI	3.3	Type 26	100425
	TFI	3.75	Type 26	100425
	TFI	4.2	Type 26	100425
	TFI	5.0	Type 26	100425
THOMMEN MEDICAL	SPI System Platform	3.5	Type 22	100404
	SPI System Platform	4.0	Type 21	100402
	SPI System Platform	4.5	Type 15	100386
	SPI System Platform	5.0	Type 15	100386
	SPI System Platform	6.0	Type 15	100386
TIDAL SPIRAL		3.75	Type 27	100431
		4.2	Type 27	100431
		5.0	Type 27	100431

Implant Company	Model	Size	SmartPeg	Art. No
TIDAL SPIRAL (Continued)		6.0	Type 27	100431
TRI DENTAL IMPLANTS	TRI-Octa 4.8	3.75	Type 61	100544
	TRI-Octa 4.8	4.1	Type 61	100544
	TRI-Octa 4.8	4.7	Type 61	100544
	TRI-Vent	3.5 platform	Type 32	100440
U-IMPL	External Hex		Type 01	100353
	Internal Hex		Type 32	100440
WARANTEC	Oneplant Implant	4.3	Type 24	100412
	Oneplant Internal Tapered		Type 06	100378
WEGO	3.4 x 9	3.4	Type 53	100484
	4.3 x 7	4.3	Type 38	100455
	4.3 x 9	4.3	Type 23	100411
WIELAND IMPLANT SYSTEM		3.5	Type 45	100472
		4.3	Type 45	100472
		5.0	Type 45	100472
YAMAHACHI DENTAL	μ -one HA Implant 2 piece	3.4	Type 38	100455
	μ -one HA Implant 2 piece	3.8	Type 38	100455
	μ -one HA Implant 2 piece	4.3	Type 27	100431
	μ -one HA Implant 2 piece	5.0	Type 18	100398
ZIMMER DENTAL	AdVent	4.5	Type 35	100446
	AdVent	5.7	Type 34	100444
	Screw-vent 3.5	3.3	Type 32	100440
	Screw-vent 3.5	3.7	Type 32	100440
	Screw-vent 4.5	4.7	Type 27	100431
	Tapered Screw-vent 3.5	3.7	Type 32	100440
	Tapered Screw-vent 3.5	4.1	Type 32	100440
	Tapered Screw-vent 4.5	4.7	Type 27	100431
	Tapered Screw-vent 5.7	6.0	Type 33	100442
	Spline	3.25	Type 30	100436
	Spline	3.75	Type 30	100436
	Spline	5.0	Type 30	100436
	SwissPlus	3.8	Type 31	100438
	SwissPlus	4.8	Type 04	100350
ZUGA MEDICAL	Bone Level	3.5	Type 63	100554
	Bone Level	5.0	Type 65	100557

Abutment System	Model	Size	SmartPeg	Art. No
COCHLEAR	Baha BA 210	5.5	Type 55	100487
	Baha BA 210	8.5	Type 55	100487
	Baha BA 300, BIA 300	6.0	Type 55	100487
	Baha BA 300, BIA 300	9.0	Type 55	100487
	Baha BA 400, BIA 400	6.0	Type 55	100487
	Baha BA 400, BIA 400	8.0	Type 55	100487
	Baha BA 400, BIA 400	10.0	Type 55	100487
	Baha BA 400, BIA 400	12.0	Type 55	100487
	Baha BA 400, BIA 400	14.0	Type 55	100487
	Vistafix Standard abutment	3.0	Type A3	100363
	Vistafix Standard abutment	4.0	Type A3	100363
	Vistafix Standard abutment	5.5	Type A3	100363
	Vistafix Standard abutment	7.0	Type A3	100363
	Vistafix VXA 300	3.5	Type A3	100363
	Vistafix VXA 300	4.5	Type A3	100363
	Vistafix VXA 300	6.0	Type A3	100363
	Vistafix VXA 300	7.5	Type A3	100363
NEOSS	Access abutment		Type 52	100483
NOBEL BIOCARE	Brånemark system MUA	NP	Type A3	100363
	Brånemark system MUA	RP	Type A3	100363
OTICON MEDICAL	Ponto Implant abutment	4.0	Type 55	100487
THOMMEN MEDICAL	SPI System Bar abutment		Type 36	100447
	VARIOMulti abutment		Type 25	100415

SmartPeg

A fresh start every time!

To avoid damage to the implant, the SmartPeg material is quite soft, which limits the lifetime of the threads. Therefore, SmartPegs are single use items.



www.osstell.com/webshop

 OSSTELL