

Osstell ISQ guides your implant treatments



Stampgatan 14, 411 01 Göteborg, Sweden
Phone: +46 31 340 82 50
info@osstell.com | www.osstell.com

Osstell Beacon Brochure
Osstell AB Date of issue: 2023-11-08 25103-02 US



OSSTELL
Beacon

Go from **Good** to **Great**

Intuitive | Easy to use | Hassle-free

The Osstell Beacon tells you when an implant is ready for loading in a matter of seconds.



PREDICTABLE OUTCOMES

Our method eliminates chance. When measuring implant stability, you will get an objective value, making implant treatments more predictable - even for patients with risk factors.




REDUCED TREATMENT TIME

Immediate loading is becoming more and more common. Patients want results fast. With Osstell®, you only wait as long as needed.



PROVEN METHOD

Backed by more than 1400 scientific studies, our method is based on RFA to determine whether or not an implant is stable enough. The result is presented as an ISQ value of 1-100. The higher the ISQ, the more stable the implant.



The Osstell ISQ Scale

The Osstell Beacon uses proprietary and evidenced-based RFA technology, backed by more than 1400 scientific studies. Visit osstell.com/scientific-database to learn more.

OsstellConnect

[OsstellConnect](#) is a cloud-based solution offering a new way for clinicians to use data to enhance their dental implant treatments. It is the world's largest implant online portal with more than a quarter of a million datapoints. OsstellConnect is currently free of charge and accessible to Osstell Beacon and Osstell IDx.

Introducing the new OsstellConnect App

Optimized for mobile and tablet. Scan your implant box for instant implant details, simplify the consultation process, access clinic data insights, and effortlessly record ISQ readings with your Beacon directly connected to your mobile device.

